

EXPERIMENT NO-6 Determination of Hardness in a sample of water:-

Background: Hard Water: Hard waters are generally considered to be those waters that require considerable amounts of soap to produce foam and that also produce scale in water pipes, heaters, boilers and other units in which the temperature of water is increased. Hard water are appropriate for human consumption similar to that as soft waters, however it produces adverse actions with soap and thus their use for cleaning purposes is unsatisfactory and thus their removal from water is required. Hardness of waters varies from place to place. In general, surface waters are softer than ground waters. Waters are commonly classified based on degree of hardness

Classification of hardness types

Hardness (mg/L)	Degree of hardness
0-75	Soft
75-100	Moderately hard
150-300	Hard
>300	Very hard

Hardness is most commonly measured by titration with an EDTA solution. A titration involves adding small amounts of a solution to a **water** sample until the sample changes color. You can titrate a sample for total **hardness** using a burette .

Total hardness is defined as the sum of the calcium and magnesium concentrations, both expressed as calcium carbonate in mg/L. When hardness (numerically) is greater than the sum of carbonate and bicarbonate alkalinity, amount of hardness equivalent to the total alkalinity is called “Carbonate hardness”. . When the hardness is numerically equal to or less than the sum of carbonate and bicarbonate alkalinity all of the hardness is carbonate hardness and there is no non carbonate hardness. The hardness may range from zero to hundreds of milligrams per litre in terms of calcium carbonate, depending on the source and treatment to which the water has been subjected.

Ethylenediamine tetra-acetic acid and its sodium salts (EDTA) form a chelated soluble complex when added to a solution of certain metal cations. If a small amount of a dye such as Eriochrome black T is added to an aqueous solution containing calcium and magnesium ions at a pH of 10 ± 0.1 , the solution will become wine red. If EDTA is then added as a titrant, the calcium and magnesium will be complexed. After sufficient EDTA has been added to complex all the

magnesium and calcium, the solution will turn from wine red to blue. This is the end point of the titration.

Apparatus

1. Burette
2. Pipette
3. Erlenmeyer flask
4. Bottle etc.

Reagents

1. Standard EDTA titrant (0.01 M)
2. Eriochrome black T indicator
3. Ammonia buffer solution

Procedure

1. Dilute 25 mL of sample (V) to about 50 mL with distilled water in an Erlenmeyer flask.
2. Add 1 mL of buffer solution.
3. Add two drops of indicator solution. The solution turns wine red in colour.
4. Add the standard EDTA titrant slowly with continuous stirring until the last reddish tinge disappears from the solution. The colour of the solution at the end point is blue under normal conditions.
5. Note down the volume of EDTA added (V_1).

$$\text{Hardness as CaCO}_3 = \frac{1000 V_1}{V} = \dots\dots\dots \text{mg / L}$$

See the link below for demonstration:-

https://www.youtube.com/watch?v=4X9r95OyOaI&ab_channel=DCVMUDCVMU